# HALOFERAX DNA PREPARATION

#### **GENOMIC DNA MINIPREP BY SPOOLING**

- Set up 5 ml liquid culture in Hv-YPC using 1-4 colonies. Grow overnight at 45°C until A650 ~0.8
- Pellet 1 ml of culture at 6,000 rpm for 8 min in round-bottomed 2 ml eppendorf. Remove supernatant
- \* Resuspend in 200 μl of <u>ST buffer</u>, scrape cells off wall of tube before pipetting (1 M NaCl, 20 mM Tris.HCl)
- Add 200 μl lysis solution to each tube and invert to mix (100 mM EDTA pH 8.0, 0.2% SDS), will go clear

#### At this point deal with each sample one at a time

- Overlay aqueous layer with 1 ml ethanol. Spool DNA at interface onto capillary tip.
  Continue spooling until liquid is homogeneous and ~clear
- Transfer spooled DNA to eppendorf with 1 ml ethanol, and swirl to wash DNA. Repeat with fresh 1 ml of ethanol. Impale pipette tip on rack to let ethanol drain from DNA while processing next sample
- Resuspend spool of DNA in 500 µl TE in eppendorf. Leave to soak in TE until all samples are processed, then scrape DNA off tip onto lip of eppendorf

#### At this point deal with all samples at the same time

- Add 50 μl <u>3 M sodium acetate</u> (pH 5.2) and 400 μl isopropanol, and invert to mix
- Pellet at Max rpm for 5-10 min. Wash pellet with 1 ml 70% ethanol and dry pellet thoroughly
- <sup>▶</sup> Resuspend DNA in 100 µl TE and incubate (shaking) at 45°C for ≥1 hour
- Leave DNA at 4°C overnight to resuspend completely. Use 5 μl per digest for Southerns

#### MAGIC BOOK

#### PLASMID MINIPREP FROM HALOFERAX VOLCANII

- Set up 10-20 ml overnight culture, depending on plasmid copy number (e.g. pHV2 vs pHV1)
- Pellet cells at 6,000 rpm for 12 minutes in round-bottom tube
- Resuspend pellet in 125 µl of <u>ST buffer</u> (1 M NaCl 20 mM Tris.HCl pH 7.5)
  - Add 25 µl of <u>500 mM EDTA</u> and invert to mix
- Add 125 µl of A1 resuspension buffer (from kit https://www.fishersci.co.uk/shop/products/ nucleospin-plasmid-columns/11932392) and mix thoroughly
- Add 250 μl of A2 lysis buffer (from kit) and mix thoroughly
  - If lysis is not complete, add more EDTA
- Add 300 µl of the neutralisation buffer A3 (from kit) and invert several times
  - Spin for 10 minutes at 12,000 rpm and transfer supernatant to fresh tube
- Add equal volume of chloroform:IAA (24:1) and vortex thoroughly
  - Spin for 1 or 2 minutes and transfer top layer to fresh tube
  - Repeat if there is a lot of white junk at the interface
- Load onto miniprep column (from kit) as usual, and follow rest of protocol as normal
  - Due to the low plasmid yield, elute DNA in small volume of elution buffer (e.g. 30 µl)

## GENOMIC DNA MAXIPREP & PHENOL-CHLOROFORM CLEANUP Day 1

Set up 5ml O/N

#### **DAY 2**

<sup>•</sup> Dilute strains 1/100 into 100ml (or more) and grow O/N until O.D. ~0.6 (or to specific O.D.)

#### DAY 3

<sup>•</sup> When the culture reaches desired O.D., spin cells down at 6,000RPM, 8 minutes, 25°C

- \* Resuspend the pellet in 5ml of <u>ST buffer</u> and transfer to a fresh 50ml Falcon tube. Add 5ml <u>lysis buffer</u>, mix and leave for ~5 minutes at RT, until clear
- <sup>•</sup> Overlay carefully with 25ml 100% ethanol, spool the DNA at the interface using a 1ml stripette, and continue spooling until homogenous and clear
- Transfer spool to a 50ml falcon containing 20ml 100% ethanol, swirl to wash. Let excess ethanol drip off until spooled DNA is semi-dry (but not completely dry!)
- <sup>•</sup> Transfer the spooled DNA to a 50ml falcon and resuspend in 15ml TE (or 5 ml, if pellet is small)
- <sup>•</sup> Add 1.5ml (or 0.5 ml, 1/10<sup>th</sup> vol) <u>3M NaAc</u> and overlay with 12ml (or 4 ml) Isopropanol and spool the DNA at the interface using a 1ml stripette, continue spooling until homogenous and clear

#### DAY 3 CONT.

- <sup>•</sup> Transfer spooled DNA to 50ml falcon with 20ml 70% ethanol, swirl to wash, repeat 70% ethanol wash. Let excess ethanol drain off until spooled DNA is semi-dry (but not completely dry)
- <sup>•</sup> Resuspend the pellet in 5 ml (or 2 ml) TE in 15 ml Falcon. Add 10 µl (or 4 µl) RNase (30mg/ml). After 1 hour at 37°C add 100 µl (or 40 µl) Proteinase K (2mg/ml in H<sub>2</sub>O). DO NOT ADD SIMULTANEOUSLY! Leave shaking gently at 37°C for 1 hour
- <sup>•</sup> Overlay DNA with equal volume of phenol:chloroform:IAA mix (25:24:1 **take care!**). Shake and leave on tube roller for 10 minutes to mix. Centrifuge for 5 minutes, 13,000 RPM, 4°C. Carefully remove the top (aqueous) layer and transfer to fresh Falcon tube, <u>avoid white interface</u>. Repeat phenol extraction
- <sup>•</sup> Overlay the DNA with equal volume chloroform:IAA mix (24:1). Shake and leave on tube roller for 10 minutes to mix. Centrifuge for 5 minutes, 13,000 RPM, 4°C. Carefully remove the top (aqueous) layer and transfer to fresh falcon tube, <u>avoiding white interface</u>. Repeat chloroform extraction <u>at least</u> once.
- <sup>•</sup> To DNA add 1/10 volume <u>3M NaAC</u> and carefully overlay with 2 volumes 100% ethanol, and spool the DNA at the interface using a gel loading tip (as in <u>genomic DNA mini prep</u>)
- Transfer spooled DNA to eppendorf with 1ml 70% ethanol, swirl to wash and repeat, then let excess ethanol drain off until spooled DNA is semi-dry (but not completely dry)
- <sup>•</sup>Resuspend in 500 µl (or 200 µl) TE, leave overnight or longer at 4°C to resuspend completely

## **ST BUFFER FOR HVO DNA PREP**

- 40 ml 5 M NaCl (to 1 M)
- 4 ml 1 M Tris.HCl pH 7.5 (to 20 mM)
- 156 ml H<sub>2</sub>0 to 200 ml
- Autoclave, store at RT

## LYSIS SOLUTION FOR HVO DNA PREP

- 40 ml 0.5 M EDTA pH 8 (to 100 mM)
- ► 4 ml 10% SDS (to 0.2%)



- 156 ml H<sub>2</sub>0 to 200 ml
- Filter sterilise (0.22  $\mu m$ ) and store at 30°C

## **0.5M EDTA (TO MAKE 1L @ PH 8.0)**

- 186.1g EDTA
- ► 800ml dH<sub>2</sub>O
- Adjust to pH 8.0 with NaOH
  - First with ~ 20g pellets added slowly followed by 5M stock solution

#### NOTE: EDTA will not go into solution until pH 8.0 is reached

#### 3M SODIUM ACETATE (@ <sub>P</sub>H 5.2)

- ▶ 61.25 g sodium acetate.3H<sub>2</sub>O
- 200 ml dH<sub>2</sub>O
- Adjust to pH 5.2 with conc. acetic acid



- ▶ make up to 250 ml with dH<sub>2</sub>O
- Filter sterilise (0.22 μm) and aliquot into 25 ml portions