


HALOFERAX TRANSFORMATION

DAY 1

- ▶ Ensure adequate quantities of transformation solutions:
 - ▶ [buffered](#) and [unbuffered](#) spheroplasting solutions
 - ▶ [spheroplast dilution solution](#)
 - ▶ [regeneration solution](#)
 - ▶ [plating solution](#)
- ▶ Need 3 selective and 1 non-selective plates per transformation
- ▶ Set up 5 ml overnight(s) in Hv-YPC (+thy), use 1–4 colonies.
 - ▶ 5 ml overnight culture normally enough for up to 3 transformations
- ▶ Leave PEG 600 at RT (or 30°C if necessary) to thaw overnight
 - ▶ Warm to 37°C if still solid


DAY 2

- ▶ When $A_{650} \approx 0.8$, pellet cells in round-bottomed 14 ml tubes at 4,400 rpm (3,000 X g) in swing-bucket centrifuge for 10 minutes at 25°C
- ▶ Remove supernatant with pipette and resuspend pellet gently in 2 ml [buffered spheroplasting solution](#). Pipette up-and-down slowly, no more than 3 times
- ▶ Transfer to 2 ml round-bottomed tube and pellet cells at 6,000 rpm (3,600 X g) in swing-bucket microfuge for 8 minutes at 25°C
- ▶ Resuspend very gently in 400 or 600 μ l of [buffered spheroplasting solution](#). Pipette up-and-down slowly, no more than 3 times, and avoid air bubbles
 - ▶ For 4 transformations use 800 μ l, i.e. 200 μ l/transformation
- ▶ For each transformation, transfer 200 μ l cells to 2 ml round-bottomed tube.
- ▶ Add 20 μ l drop of 0.5 M EDTA pH 8.0 on side of tube, then invert to mix. Leave at RT for **10 minutes** to form spheroplasts.



DAY 2 CONT.

- ▶ Meanwhile set up DNA samples in 30 μl total:
 - ▶ 10 μl of dam⁻ DNA (~1–2 μg , to 10 μl with dH₂O) or 10 μl of dH₂O (control)
 - ▶ 15 μl of [unbuffered spheroplasting solution](#)
 - ▶ 5 μl of 0.5 M EDTA pH 8.0
- ▶ After 10 minutes add DNA in same manner as EDTA.
Leave at RT for **5 minutes**
- ▶ Meanwhile prepare 60% PEG 600 solution. For 3 transformation reactions use:
 - ▶ 480 μl of PEG 600 (viscous, pipette very slowly for accuracy)
 - ▶ 320 μl of [unbuffered spheroplasting solution](#).
- ▶ After 5 minutes add 250 μl (equal volume) of 60% PEG 600 to each transformation. Add in same manner as EDTA, but shake tube horizontally to mix gently. Leave at RT for **30 minutes**
 - ▶ 1 hour for *H. mediterranei*



DAY 2 CONT.

- ▶ Add 1.5 ml [spheroplast dilution solution](#), invert to mix and leave at RT for **2 minutes**
- ▶ Pellet cells in swing-bucket microfuge at 6,000 rpm (3,600 X g) for 8 minutes at 25°C, and remove supernatant with pipette
- ▶ Next to flame and using filter tip, add 1 ml [regeneration solution](#) (+ 60 µg/ml thymidine, if required)
- ▶ Use blue filter tip to scrape cells off tube wall, suck up pellet (but avoid resuspending pellet), then transfer to 4 ml sterile tube
- ▶ Put tube in rotator at 45°C for **4 hours or longer**
 - ▶ Not necessary to leave undisturbed for 1 hour first

DAY 2 CONT.

- ▶ Transfer cells to 2 ml round-bottomed tube, pellet in swing-bucket microfuge at 6,000 rpm (3,600 X g) for 8 minutes at 25°C
- ▶ Remove supernatant and resuspend gently in 1 ml [plating solution](#)
 - ▶ This step is not necessary if NovR or MevR selection is used.
For Thy+ selection, [regeneration solution](#) can be used
- ▶ Dilute in [plating solution](#) with following dilutions and use glass spreader to plate 100 µl on Hv-Ca (or Hv-YPC or Hv-Min) plates:
 - ▶ Use [regeneration solution](#) for NovR or MevR selection

	Selective	Non-selective (optional)
+DNA	Use 2-3 dilutions depending on transforming DNA 10 ⁻¹ , 10 ⁻³ , or 10 ⁰ , 10 ⁻² , 10 ⁻⁴ (HVO, HME)	10 ⁻⁶
No DNA control (optional)	Use 2 dilutions depending on selection 10 ⁰ , 10 ⁻¹	10 ⁻⁶

- ▶ Leave plates at 45°C for at least 5 days

BUFFERED SPHEROPLASTING SOLUTION

- ▶ 14.61g NaCl (1 M)
- ▶ 0.5g KCl (27mM)
- ▶ 12.5 ml 1M Tris.HCl pH 8.5 (50 mM)
- ▶ 37.5g Sucrose (15%)
- ▶ dH₂O to 250ml
- ▶ Filter sterilise

UNBUFFERED SPHEROPLASTING SOLUTION

- ▶ 5.84g NaCl (1 M)
- ▶ 0.2g KCl (27 mM)
- ▶ 15g Sucrose (15%)
- ▶ dH₂O to 100ml
- ▶ adjust to pH 7.5 (~10µl 1M NaOH)
- ▶ Filter sterilise

SPHEROPLAST DILUTION SOLUTION

- ▶ 76ml 30% SW (23%)
- ▶ 15g Sucrose (15%)
- ▶ 0.75ml 0.5M CaCl₂ (3.75 mM)
- ▶ dH₂O to 100ml
- ▶ Filter sterilise

REGENERATION SOLUTION

- ▶ 150ml 30% SW (18%)
- ▶ 25ml 10X YPC (1X)
- ▶ 37.5g Sucrose (15%)
- ▶ 1.5ml 0.5M CaCl₂ (3 mM)
- ▶ dH₂O to 250ml
- ▶ Filter sterilise

PLATING SOLUTION

- ▶ 150ml 30% SW (18%)
- ▶ 37.5g Sucrose (15%)
- ▶ 1.5ml 0.5M CaCl₂ (3 mM)
- ▶ dH₂O to 250ml
- ▶ Filter sterilise