

HALOFERAX MATING

HALOFERAX MATING

DAY 1



- Set up 5ml O/Ns of strains to be mated in Hv-broth (YPC or Ca, +Thy if needed)

DAY 2

- Dilute O/Ns in fresh broth, use 10 μ l (1/500) for 5ml O/N
- Prepare [Millipore Swinnex 25 filter units](#) with [Millipore 0.45 \$\mu\$ m type MCE 25 mm filters](#), wrap in aluminium foil and autoclave

DAY 3


- When OD = 0.8 ($\sim 4 \times 10^8$ cells) for both strains, mix equal volumes (1 ml each) of O/Ns
 - Also use 2 ml of unmixed parental O/Ns as negative controls
 - For different OD values, see this [table](#)
- Use 5 ml syringe to gently filter cell suspension through the sterile Swinnex filter unit
 - Repeat with 2 ml of unmixed parental O/Ns as negative controls
- Disassemble filter unit and use sterile forceps to remove filter disc. Place filter cell-side uppermost on Hv-YPC (+Thy) plate)
- Incubate overnight at 45°C



HALOFERAX MATING CONTINUED...


DAY 4

- Use sterile forceps to remove filter disc, transfer to 5 ml eppendorf with 1 ml 18% SW
 - Leave in rotator at 45°C ~1 hour to gently resuspend cells
- Dilute cells in 18% SW, plate 100 µl at 10⁰, 10⁻¹ and 10⁻² on selective (e.g. Hv-Ca) plates
 - Optional: spot 20 µl at 10⁻⁴, 10⁻⁶ and 10⁻⁸ on Hv-YPC for viable count
 - Incubate at 45°C for 5-6 days; expect recombinants at frequency of 10⁻⁴



DAY 10 ONWARDS


- Set up 5ml O/Ns of heterodiploid cells in non-selective broth (e.g. Hv-YPC, +Thy)
 - Optional: use [colony PCR](#) to verify heterodiploid cells
 - e.g. [PCR for polB](#), differs between Hvol and Hmed (intein in former), expect both bands in diploid



- Next day, dilute cells in 18% SW, plate 100 µl at 10⁻¹, 10⁻², 10⁻³ and 10⁻⁴ on selective plates
 - e.g. Hv-Ca+5FOA, incubate at 45°C for 5-6 days


- Recombinant haploids will have novel combination of selectable markers

- Optional: use [colony PCR](#) to verify haploid recombinants, e.g. *polB*, expect one band only



CELLS PER ML @ OD650 VALUE

OD650	Cells/ml
0.022	9.00E+05
0.05	3.00E+06
0.09	1.20E+07
0.15	2.40E+07
0.25	8.00E+07
0.39	1.40E+08
0.86	3.50E+08
0.1	2.50E+07
0.2	5.00E+07
0.3	7.50E+07
0.4	1.00E+08
0.5	1.75E+08
0.6	2.50E+08
0.7	3.25E+08
0.9	4.00E+08



COLONY PCR: ONETAQ

DNA prep:

- ▶ Touch DNA with yellow tip (the less cells the better)
- ▶ Pipette cells up and down in 100µl SDW
- ▶ Boil at 100°C for 10mins
- ▶ Leave on ice for 10mins

PCR mix per Reaction	
DNA	1
Forward primer	0.4
Reverse primer	0.4
dNTPs	4
GC buffer (5X)	4
OneTaq	0.1
SDW	10.1
20µl	

PCR cycles	Temperature /°C	Time	x 30
Initial D	94	30s	
D	94	30s	
A	x	30s	
E	68	y (1min per kb)	
Final E	68	5mins	

Can use [Q5 hotstart polymerase](#) instead, but costs more

Q5 PCR

PCR mix per Reaction	
DNA*	1
Forward primer (50 µM)	1
Reverse primer (50 µM)	1
dNTPs (1 mM)	10
GC enhancer (5X)	10
Q5 buffer (5X)	10
Q5 hotstart polymerase	0.5
SDW	16.5
50µl	

- ▶ DNA* Use ~50 ng plasmid DNA; if using genomic DNA, dilute 1/10 in dH₂O
- ▶ x** = [annealing temperature](#)
= 81.5 + 0.41%GC - 600/Length + 16.6Log₁₀[Na⁺]
Assuming Na⁺ concentration = 50 mM
Subtract 1°C for each 1% mismatch
- ▶ y[†] = extension time, 20 sec/kb for plasmid DNA, 30 sec/kb for genomic DNA,

PCR cycles	Temperature /°C	Time	
Initial Denature	98	30s	
Denature	98	10s	x 30
Anneal	x**	20s	
Extension	72	y [†] (secs per kb)	
Final Extension	72	10mins	

TOUCHDOWN PCR

TD PCR cycles	Temperature /°C	Time	
Initial Denature	98	30s	
Denature	98	10s	x 10
Anneal	Low - High temp**	20s	
Extension	72	y [†] (secs per kb)	
Denature	98	10s	x 20
Anneal	High temp**	20s	
Extension	72	y [†] (secs per kb)	
Final Extension	72	10mins	