


HALOFERAX MEDIA

MARKERS VS MEDIA



Marker	HV-YPC	HV-CA	HV-MIN
<i>pyrE2</i>		+ura	+ura
<i>trpA</i>		+trp	+trp
<i>leuB</i>			+leu
<i>hdrB</i>	+thy	+thy +hypox	+thy +hypox +met +gly +panto

30% SALT WATER (30% SW)

- ▶ 4L dH₂O, microwave until hot and stir on hot plate
- ▶ 1200g NaCl
- ▶ 150g MgCl₂·6H₂O
- ▶ 175g MgSO₄·7H₂O
- ▶ 35g KCl
- ▶ 100ml 1M Tris.HCl pH7.5
- ▶ Make up to 5L with dH₂O
 - ▶ To make 20% SW, use 2L 30% SW and 1L dH₂O

10X YPC (FOR 5 BOTTLES)

- ▶ 130ml dH₂O
- ▶ 8.5g yeast extract
- ▶ 1.7g peptone (Oxoid)
- ▶ 1.7g casamino acids
- ▶ 3ml 1M KOH (add slowly)
- ▶ Make up 170ml with dH₂O

10X YPC (FOR 10 BOTTLES)

- ▶ 250ml dH₂O
- ▶ 17g yeast extract
- ▶ 3.4g peptone (Oxoid)
- ▶ 3.4g casamino acids
- ▶ 6ml 1M KOH (add slowly)
- ▶ Make up 340ml with dH₂O

HV-YPC BROTH

- ▶ Into each 500ml Duran bottle (media grade) measure:
- ▶ 300ml [20% SW](#)
 - ▶ Or 200ml [30% SW](#) and 100ml dH₂O
- ▶ 33ml [10X YPC](#)
- ▶ Autoclave for 1 minutes at 121°C (program 1)
- ▶ When cool, add 2ml [0.5M CaCl₂](#) and store in dark cupboard

HV-YPC BROTH – LARGE VOLUMES

	Molarity	1 litre	2 litres	4 litres	8 litres
dH ₂ O		800 ml	1.600 ml	3.200 ml	6.400 ml
NaCl	2,4641 M	144 g	288 g	576 g	1.152 g
MgSO ₄ .7H ₂ O	0,0852 M	21 g	42 g	84 g	168 g
MgCl ₂ .6H ₂ O	0,0885 M	18 g	36 g	72 g	144 g
KCl	0,0563 M	4 g	8 g	17 g	34 g
1M Tris pH7.5	0,012 M	12 ml	24 ml	48 ml	96 ml
Yeast extract	0.5% (wt/vol)	5 g	10 g	20 g	40 g
Peptone	0.1% (wt/vol)	1 g	2 g	4 g	8 g
Casamino acids	0.1% (wt/vol)	1 g	2 g	4 g	8 g
Add dH ₂ O to		1.000 ml	2.000 ml	4.000 ml	8.000 ml
Sterilise in large autoclave on sugar run and leave to cool overnight					
0.5 M CaCl ₂	0,003 M	6 ml	12 ml	24 ml	48 ml

HV-YPC AGAR

- ▶ Into each 500ml Duran bottle (media grade) measure:
 - ▶ 300ml [20% SW](#)
 - ▶ 5g agar
- ▶ Microwave, for 10 bottles: 20 min on High (100%), 20 min on Medium/High (70%), 2x 10 min on Medium (50%)
- ▶ Add 33ml [10X YPC](#) and large stir bar
- ▶ Autoclave immediately for 1 minutes at 121°C (program 1)
- ▶ Leave in water bath at 60°C to cool

HV-YPC AGAR PLATES

- ▶ Mark plates on side of lid with thick **BLUE** stripe
- ▶ Use thin stripe for optional additives and write on plate:
 - ▶ 3.4ml **thymidine** (1/100 of 4mg/ml in dH₂O, to 40µg/ml), **red**
 - ▶ 67µl **novobiocin** (1/500 of 1mg/ml in dH₂O, to 0.2µg/ml), **black**
 - ▶ 67µl **mevinolin** (1/5000 of 10mg/ml in ethanol, to 2µg/ml), **purple**
- ▶ When cool, add 2ml 0.5M CaCl₂ plus additives while stirring
- ▶ Pour 9 thick plates of ~37ml per bottle
- ▶ Seal 18 plates per bag, in **YPC** box in cold room, update board

10X CA (FOR 5 BOTTLES)

- ▶ 130ml dH₂O
- ▶ 8.5g casamino acids
- ▶ 3ml 1M KOH (add slowly)
- ▶ Make up 170ml with dH₂O

10X CA (FOR 10 BOTTLES)

- ▶ 250ml dH₂O
- ▶ 17g casamino acids
- ▶ 6ml 1M KOH (add slowly)
- ▶ Make up 340ml with dH₂O

HV-CA BROTH

- ▶ Into each 500ml Duran bottle (media grade) measure:
 - ▶ 200ml [30% SW](#)
 - ▶ 10ml 1M Tris.HCl [pH7.0](#)
 - ▶ 77ml dH₂O
 - ▶ 33ml [10X Ca](#)
- ▶ Autoclave for 1 minutes at 121°C (program 1)
- ▶ When cool, add:
 - ▶ 8.5ml [Hv-Min Carbon source](#)
 - ▶ 4ml [Hv-Min salts](#)
 - ▶ 666µl [0.5M KPO₄ buffer](#) (pH7.0)
 - ▶ 435µl [thiamine&biotin](#)
- ▶ Store in dark cupboard

HV-CA AGAR

- ▶ Into each 500ml Duran bottle (media grade) measure:
 - ▶ 300ml [20% SW](#)
 - ▶ 5g agar
- ▶ Microwave, for 10 bottles: 20 min on High (100%), 20 min on Medium/High (70%), 2x 10 min on Medium (50%)
- ▶ Add 33ml [10X Ca](#) and large stir bar
- ▶ Autoclave immediately for 1 minutes at 121°C (program 1)
- ▶ Leave in water bath at 60°C to cool

HV-**Ca** AGAR PLATES

- ▶ Mark plates on side of lid with thick **RED** stripe
- ▶ Use thin stripe for optional additives and write on plate:
 - ▶ 3.4ml **thymidine&hypoxanthine** (1/100 of 4mg/ml ea. in 50 mM NaOH to 40µg/ml), **red**
 - ▶ 1.7ml **tryptophan** (1/200 of 10mg/ml in dH₂O to 50 µg/ml), **brown**
 - ▶ 340µl **uracil** (1/1000 of 50mg/ml in DMSO to 50µg/ml), **pink**
 - ▶ 340µl **5-FOA** (1/1000 of 50mg/ml in DMSO to 50µg/ml) +67µl uracil (to 10µg/ml), **2x pink**
- ▶ When cool, add while stirring gently:
 - ▶ 2.8ml Hv-Ca salts
 - ▶ 666µl 0.5M KPO₄ buffer (pH7.0)
 - ▶ Plus additives
- ▶ Pour 9 thick plates of ~37ml per bottle
- ▶ Seal 18 plates per bag, store in **Ca** box in cold room, update board

HV-MIN BROTH – UREA-BASED

- ▶ Into each 500ml Duran bottle (media grade) measure:
 - ▶ 200ml [30% SW](#)
 - ▶ 10ml 1M Tris.HCl [pH7.0](#)
 - ▶ 77ml dH₂O
- ▶ Autoclave for 1 minutes at 121°C (program 1)
- ▶ When cool, add:
 - ▶ 8.5ml [Hv-Min carbon source](#)
 - ▶ 4ml [Hv-Min salts \(urea\)](#)
 - ▶ 666µl 0.5M [KPO4 buffer](#) (pH7.0)
 - ▶ 435µl [thiamine&biotin](#)
- ▶ Store in dark cupboard

HV-“MIN” BROTH – AMINO ACID-BASED

- ▶ Into each 500ml Duran bottle (media grade) measure:
 - ▶ 200ml [30% SW](#)
 - ▶ 10ml 1M Tris.HCl [pH7.0](#)
 - ▶ 77ml dH₂O
- ▶ Autoclave for 1 minutes at 121°C (program 1)
- ▶ When cool, add:
 - ▶ 35ml [Hv-Min amino acid mix](#)
 - ▶ 8.5ml [Hv-Min carbon source](#)
 - ▶ 2ml [0.5M CaCl₂](#)
 - ▶ 333µl [trace elements](#)
 - ▶ 666µl [0.5M KPO₄ buffer](#) (pH7.0)
 - ▶ 435µl [thiamine&biotin](#)
- ▶ Store in dark cupboard

HV-MIN AGAR

- ▶ Into each 500ml Duran bottle (media grade) measure:
 - ▶ 300ml [20% SW](#)
 - ▶ 10ml dH₂O
 - ▶ 10ml 1M Tris.HCl [pH7.0](#)
 - ▶ 5g agar
- ▶ Microwave, for 10 bottles: 20 min on High (100%), 20 min on Medium/High (70%), 2x 10 min on Medium (50%)
- ▶ Add large stir bar, autoclave immediately for 1 minutes at 121°C
- ▶ Leave in water bath at 60°C to cool

HV-MIN AGAR PLATES

- ▶ Mark plates on side of lid with thick **GREEN** stripe
- ▶ Use thin stripe for optional additives and write on plate:
 - ▶ 3.4ml **thymidine&hypoxanthine** (1/100 of 4mg/ml each in 50 mM NaOH to 40µg/ml), **red**
 - ▶ Plus 1.7ml methionine & glycine & pantothenic acid (1/200 of 10 mg/ml each in dH₂O, to 50µg/ml)
 - ▶ 1.7ml **tryptophan** (1/200 of 10mg/ml in dH₂O to 50 µg/ml), **brown**
 - ▶ 340µl **uracil** (1/1000 of 50mg/ml in DMSO to 50µg/ml), **pink**
 - ▶ 1.7ml **leucine** (1/200 of 10mg/ml in dH₂O to 50 µg/ml), **green**
 - ▶ 2.2ml **histidine** (1/150 of 10mg/ml in dH₂O to 65 µg/ml), **blue**
 - ▶ 1.7ml **methionine** (1/200 of 10mg/ml in dH₂O to 50 µg/ml), **orange**
- ▶ When cool, add while stirring gently:
 - ▶ 8.5ml Hv-Min Carbon Source
 - ▶ 4ml Hv-Min salts
 - ▶ 666µl 0.5M KPO₄ buffer (pH7.0)
 - ▶ 435µl thiamine&biotin
 - ▶ Plus additives
- ▶ Pour 9 thick plates of ~37ml per bottle, seal 18 plates per bag, store in **Min** box in cold room

80% GLYCEROL 6% SW

- ▶ 80ml glycerol
- ▶ 20ml [30% SW](#)
- ▶ Autoclave for 1 minutes at 121°C (program 1)
- ▶ When cool add 200µl [0.5M CaCl₂](#)

18% SW

- ▶ 200ml [30% SW](#)
- ▶ 133ml dH₂O
- ▶ Autoclave for 1 minutes at 121°C (program 1)
- ▶ When cool add 2ml [0.5M CaCl₂](#)

HV-CA SALTS

- ▶ 13ml [thiamine&biotin](#)
- ▶ 60ml 0.5M CaCl₂
- ▶ 10ml [trace elements](#)
- ▶ Store at 4°C

0.5M CaCl₂

- ▶ 800ml dH₂O
- ▶ 73.5g CaCl₂
- ▶ Make up 1L with dH₂O
- ▶ Filter sterilise

HV-MIN SALTS

- ▶ 30ml 1M NH₄Cl (sterile)
- ▶ 36ml 0.5M CaCl₂
- ▶ 6ml [trace elements](#)
- ▶ Store at 4°C

HV-MIN SALTS – UREA

- ▶ 30ml 1M urea (sterile)
- ▶ 36ml 0.5M CaCl₂
- ▶ 6ml [trace elements](#)
- ▶ Store at 4°C

0.5M KPO₄ BUFFER pH7.0

- ▶ 61.5ml 1M K₂HPO₄
- ▶ 38.5ml 1M KH₂PO₄
- ▶ Check pH = 7.0
- ▶ Store at 4°C

THIAMINE&BIOTIN

- ▶ 12ml thiamine
(1mg/ml in dH₂O)
- ▶ 7.5ml biotin
(0.2mg/ml in dH₂O)
- ▶ Store at 4°C

TRACE ELEMENTS

- ▶ 250ml dH₂O
 - ▶ Add a few drops cHCl
- ▶ 90mg MnCl₂·4H₂O
- ▶ 110mg ZnSO₄·7H₂O
- ▶ 575mg FeSO₄·7H₂O
- ▶ 12.5mg CuSO₄·5H₂O
- ▶ Filter sterilise
- ▶ Store at 4°C

HV-MIN CARBON SOURCE

- ▶ 150ml dH₂O
- ▶ 41.7ml 60% DL-lactic acid Na⁺salt
- ▶ 33.3g succinic acid disodium salt·6H₂O
- ▶ 6.3ml 80% glycerol
- ▶ Check pH = 7.0, make up 250 ml with dH₂O
- ▶ Filter sterilise, store at 4°C

AMINO ACID MIX (50mM)

- ▶ Weigh out amino acids
 - ▶ See table for proportions
- ▶ Blend powders in grinder
 - ▶ Pulse to avoid overheating
- ▶ Store in dark cupboard
- ▶ Dissolve 654.14mg mix in 100ml dH₂O for 50mM stock
- ▶ Filter sterilise, store at 4°C

AMINO ACID COCKTAIL

Name			MW	50 ml stock (50 mM)	100 ml Cocktail (50 mM)	100 ml Mix (50 mM)	1 litre Mix (50 mM)	5 litre Mix (50 mM)
Alanine	Ala	A	89,09	222.73 mg	10.42 ml	46.42 mg	464.2 mg	2.321 g
Cysteine	Cys	C	121,16	302.90 mg	0.69 ml	4.18 mg	41.8 mg	0.209 g
Aspartic Acid	Asp	D	133,1	332.75 mg	9.47 ml	63.02 mg	630.2 mg	3.151 g
Glutamic Acid	Glu	E	185,2	463.00 mg	9.38 ml	86.86 mg	868.6 mg	4.343 g
Phenylalanine	Phe	F	165,19	412.98 mg	3.13 ml	25.85 mg	258.5 mg	1.293 g
Glycine	Gly	G	75,07	187.68 mg	8.26 ml	31.00 mg	310.0 mg	1.550 g
Histidine	His	H	155,15	387.88 mg	2.08 ml	16.14 mg	161.4 mg	0.807 g
Isoleucine	Ile	I	131,17	327.93 mg	3.96 ml	25.97 mg	259.7 mg	1.299 g
Lysine	Lys	K	146,19	365.48 mg	2.47 ml	18.05 mg	180.5 mg	0.903 g
Leucine	Leu	L	131,17	327.93 mg	8.19 ml	53.71 mg	537.1 mg	2.686 g
Methionine	Met	M	149,21	373.03 mg	1.83 ml	13.65 mg	136.5 mg	0.683 g
Asparagine	Asn	N	132,12	330.30 mg	2.65 ml	17.51 mg	175.1 mg	0.875 g
Proline	Pro	P	115,13	287.83 mg	4.42 ml	25.44 mg	254.4 mg	1.272 g
Glutamine	Gln	Q	146,14	365.35 mg	2.61 ml	19.07 mg	190.7 mg	0.954 g
Arginine	Arg	R	174,2	435.50 mg	6.65 ml	57.92 mg	579.2 mg	2.896 g
Serine	Ser	S	105,09	262.73 mg	5.57 ml	29.27 mg	292.7 mg	1.463 g
Threonine	Thr	T	119,12	297.80 mg	6.04 ml	35.97 mg	359.7 mg	1.799 g
Valine	Val	V	117,15	292.88 mg	8.54 ml	50.02 mg	500.2 mg	2.501 g
Tryptophan	Trp	W	204,23	510.58 mg	0.95 ml	9.70 mg	97.0 mg	0.485 g
Tyrosine	Tyr	Y	181,19	452.98 mg	2.69 ml	24.37 mg	243.7 mg	1.219 g
					100.00 ml	654.14 mg	6541.4 mg	32.71 g

STOCK SOLUTIONS

0.5M EDTA (TO MAKE 1L @ pH 8.0)

- ▶ 186.1g EDTA
- ▶ 800ml dH₂O
- ▶ Adjust to pH 8.0 with NaOH
 - ▶ First with ~ 20g pellets added slowly followed by 5M stock solution



NOTE: EDTA will not go into solution until pH 8.0 is reached

10X TBE (TO MAKE 1L)

- ▶ 108 g Tris base
- ▶ 55 g Boric acid
- ▶ 40 ml 0.5M EDTA pH 8.0
- ▶ make up to 1 L with dH₂O
- ▶ Autoclave

50X TAE (TO MAKE 1L)

- ▶ 242 g Tris Base
- ▶ 57.1 ml conc. acetic acid
- ▶ 100 ml 0.5M EDTA pH 8.0
- ▶ make up to 1 L with dH₂O



5X DNA GEL LOADING DYE


- ▶ 2.5 ml 1 M Tris-HCl pH 7.5
- ▶ 10 ml 0.5 M EDTA
- ▶ 7.5 g Ficoll
- ▶ 0.125 g Bromophenol blue
- ▶ 0.125 g Xylene cyanol
- ▶ make up to 50 ml with SDW, store at RT

DNA LADDER

- ▶ 200 μ l 1kb ladder (500 μ g/ml, without dye)
- ▶ 100 μ l 10X buffer 4 or CutSmart
- ▶ 200 μ l 5X loading dye
- ▶ 500 μ l SDW
- ▶ Store at 4°C

3M SODIUM ACETATE (@ pH 5.2)

- ▶ 61.25 g sodium acetate.3H₂O
- ▶ 200 ml dH₂O
- ▶ Adjust to pH 5.2 with conc. acetic acid
- ▶ make up to 250 ml with dH₂O
- ▶ Filter sterilise (0.22 µm) and aliquot into 25 ml portions



1MM dNTPS

- ▶ 20 μ l of each dNTP (100 mM stock)
- ▶ 1920 μ l dH₂O
- ▶ Aliquot into 100 μ l portions, store at -20°C

4X SDS-PAGE LOADING DYE

- ▶ 0.5 ml 1 M Tris pH 6.8
- ▶ 7.25 ml dH₂O
- ▶ 1 ml 20% SDS  
- ▶ 10 mg bromophenol blue
- ▶ 1.25 ml 80% glycerol

20% SDS

- ▶ 200 g SDS  
- ▶ 900 ml dH₂O
- ▶ Heat and stir SLOWLY until dissolved
- ▶ make up to 1 L with dH₂O

NOTE: do not shake as SDS will foam up