

PULSED FIELD GEL

PFG – SOLUTIONS FOR DNA PREP IN AGAROSE PLUGS

Lysis Solution for Plugs

- ▶ 1.5 ml 3 M Tris.HCl pH 8.8 (~20 mM)
250 ml 0.5 M EDTA pH 8 (~500 mM)
2.5 g N-lauroylsarcosine (1%)
Total ~250 ml
No need to sterilise (already sterile)
- ▶ If required, add immediately before:
Add RNase (30 mg/ml) to 10 µg/ml
Add proteinase K (powder) to 1 mg/ml

DAY 1

- ▶ Set up overnight culture in 5 ml Hv-YPC

Wash Solution for Plugs

- ▶ 6.25 ml 1 M Tris.HCl pH 7.5 (25 mM)
50 ml 0.5 M EDTA pH 8 (100 mM)
dH₂O to 250 ml
Autoclave
- ▶ If required:
Add 1 mM PMSF (100 mM in ethanol)



Wash Solution + Glycerol for Plug Storage

- ▶ 2.5 ml 1 M Tris.HCl pH 7.5 (25 mM)
20 ml 0.5 M EDTA pH 8 (100 mM)
50 ml glycerol (50%)
Autoclave
dH₂O to 100 ml

DAY 2

- Prepare 1% low-melt agarose ([SeaPlaque GTG](#)) in 18% SW:
 - Add 0.1 g low-melt agarose to 4 ml dH₂O, boil to melt
 - Add 6 ml hot (>80°C) 30% SW (only required for *H. volcanii* – use dH₂O for e.g. *E. coli*)
 - Mix thoroughly, boil once more and leave at 42°C (for 30 minutes) to cool
- Pellet 1 ml cells at 6000 rpm for 10 minutes, use 2 ml round-bottom tube
Use 2 ml for slow-growing strains (e.g. *ΔradA*)
- Resuspend cell pellet gently in 100 μl of 18% SW
 - Transfer 80 μl of cells to fresh tube, leave at 42°C for 5 minutes
- Add 320 μl of 1% low-melt agarose at 42°C, and mix gently by pipetting.
Immediately pipette 80 μl into each plug mould, avoiding air bubbles
 - Leave in fridge for 10 minutes, normally get 3-4 plugs per sample
- Transfer agarose plug to 2 ml tube containing 2 ml lysis solution +proteinase K (1 mg/ml)
 - Incubate shaking at 45°C for 3-4 hours (low speed), can put 2 plugs in each tube
 - Clean plug moulds: rinse in hot water to remove agarose, rinse in IMS, wash in dH₂O, leave to air dry
- Remove lysis solution, add 2 ml fresh lysis solution +proteinase K (1 mg/ml) +RNase (10 μg/ml)
 - Incubate shaking at 45°C overnight

DAY 3 – OPTION 1 – GAMMA RADIATION

- ▶ Remove lysis solution with pipette, replace with 2 ml wash solution. Incubate shaking at 37°C for 30 minutes
- ▶ Transfer plugs into fresh 2 ml tubes containing 0.5 ml wash solution
- ▶ Gamma irradiate plugs for required exposure time (optional)
 - ▶ [Beverley \(1989\)](#) estimates rate of 7×10^{-6} DSB / kb / Gy
 - ▶ ^{137}Cs source, 30.17 yr half-life: 9.2333 Gy/min (Jan 1990) = 4.13175 Gy/min (Jan 2025)
 - ▶ For *H. volcanii*, normally ~5-10 mins for chromosome, ~10-20 mins for pHV3 or pHV1
- ▶ Transfer plugs to fresh 2 ml tubes with 2 ml wash solution +PMSF (1 mM). Incubate with shaking at 37°C for 1 hour.
 - ▶ Repeat wash twice more, once with PMSF and once without
 - ▶ 100 mM stock solution of PMSF in 100% ethanol (8.7 mg in 500 μ l), needed to inactivate proteinase K



DAY 3 – OPTION 2 – RESTRICTION DIGEST

- ▶ Remove lysis solution with pipette and replace with 2 ml wash solution.
Incubate shaking at 37°C for 30 minutes
- ▶ Transfer agarose plugs to fresh 2 ml tubes containing 2 ml wash solution + PMSF (1 mM).
Incubate shaking at 37°C for 1 hour
 - ▶ *100 mM stock solution of PMSF in 100% ethanol (8.7 mg in 500µl), to inactivate proteinase K*
- ▶ Repeat wash step one more time with PMSF and one without PMSF, then:
 - ▶ Wash plugs 2X in 2 ml wash solution diluted 1/10 for 30 minutes at 37°C
 - ▶ Wash plugs 1X in 2 ml wash solution diluted 1/100 for 30 minutes at 37°C
 - ▶ Wash plugs 1X in 2 ml wash solution diluted 1/1000 for 30 minutes at 37°C
- ▶ Transfer plugs to fresh 1.5 ml tubes containing 1 ml 1X restriction buffer
Incubate shaking at 37°C for 1 hour
- ▶ Remove restriction buffer and replace with 250 µl fresh 1X restriction buffer
- ▶ Add 5 µl restriction enzyme (50 U, adjust accordingly), incubate shaking at 37°C overnight



DAY 3 – POURING AND LOADING GEL

- ▶ For temporary storage of plugs, transfer into 1 ml fresh wash solution, leave at 4°C
 - ▶ For long term storage, equilibrate plugs in 1 ml wash solution +50% glycerol, then store at -20°C
- ▶ Prepare 1.2% PFG gel, ([SeaKem Gold agarose](#)) ★
 - ▶ For small gel melt 1.26 g PFG agarose in 105 ml 0.5X TBE. Pour 100 ml gel
 - ▶ For large gel melt 1.98 g PFG agarose in 165 ml 0.5X TBE. Pour 160 ml gel
 - ▶ Allow to set for >1 hour (keep ~5 ml molten agarose for later)
- ▶ Wash plugs twice in 2 ml 0.5X TBE buffer at 37°C for 30 mins before loading
- ▶ Set up PFG apparatus ([BioRad CHEF Mapper](#)) with 2.2 litres of 0.5X TBE
 - ▶ Leave circulating at 14°C, for 30 minutes to cool
- ▶ Insert plug into each well (between 1/3 and whole plug, depending on well size)
 - ▶ Seal plugs in place with molten agarose (from earlier) – use as little as possible
 - ▶ Use thin slice of [lambda ladder](#) as marker – also seal in place with molten agarose
- ★ *Alternatively, arrange plugs on teeth of PFG comb (horizontal), seal in place with molten agarose along top edge, place comb in PFG mould, and pour agarose gel around comb*
 - ▶ *Remove comb carefully (watch plugs), insert slice of lambda ladder in well, fill empty wells with agarose*

DAY 3 – POURING AND LOADING GEL – CONT...

- ▶ Place gel into frame in PFG apparatus, leave to equilibrate for ~30 minutes at 14°C with circulating buffer
 - ▶ Watch out for air bubbles in circulation tubing
- ▶ Run 1.2% PFG at 14°C with following conditions:
 - ▶ Running time of 20 hrs 46 mins (vary as appropriate)
 - ▶ Two state, 120° included angle
 - ▶ Gradient 6 V/cm (0 ramping)
 - ▶ Initial switch time 0.64 sec
 - ▶ Final switch time 1 min 13.22 sec



DAY 4/5

- ▶ Stain gel with ethidium bromide and visualise
- ▶ Proceed to Southern blot protocol if required
 - ▶ Acid nick, denature and blot for 1.5x longer than usual (due to higher agarose%)

