

# RADIATION

### DAY 1


- ▶ Make pre-hyb and hyb solutions
- ▶ **Prehyb**
- ▶ Sign into workstation. Boil 800µl salmon sperm DNA (green lid) for 5 mins at 100°C, then to pre-hyb solution
- ▶ Pour pre-hyb solution into warmed hyb tube and use to 'stick' membranes onto sides of hyb tube.
- ▶ Incubate at 65°C for >3 hours (shorter is OK for colony hyb)

## DAY 1 CONTINUED



- ▶ **Hyb**
- ▶ Label all tubes if using more than one probe
  - ▶ Add 30µl TE to a 1.5ml eppendorf – keep on ice  
**For Colony Hyb** – Add 3µl probe DNA and 12µl SDW to a screw cap tube, keep on ice  
**For Southern Blot** – Add 3µl probe DNA, 3µl diluted 1kb ladder (1ng/µl) and 8µl SDW to a screw cap tube, keep on ice
  - ▶ Get 450µl salmon sperm DNA (red lid) out of freezer and keep on ice
  - ▶ Boil probe DNA for 5mins at 100°C then keep on ice – pulse spin.
  - ▶ Add 4µl [HiPrime](#) and 2µl  $^{32}\text{P}$ -dCTP (1µl for colony hyb) to DNA, mix gently and incubate at 37°C for 15-20mins.


## DAY 1 CONTINUED



### ▸ Hyb continued

- Get [gel column](#) from fridge – ‘flick’ tube to get packing buffer out of the lid. Spin out packing buffer from column: snap bottom off, take the lid off, put in collection tube and spin for 4 mins at 3600 rpm
- When incubation is finished, add the 30µl TE to probe before loading onto column (column minimum volume is 50µl).
- Spin column for 4 mins at 3600 rpm (with lids on) to purify probe away from unincorporated dCTPs – probe is the flow through!
- Add probe to 450µl salmon sperm DNA and boil for 5 mins at 100°C, then put on ice
- Pour pre-hyb solution down the sink (with lots of water) out of small hole in hyb tube
- Use syringe barrel to add hyb solution to hyb tube. Then use p1000 pipette to add probe vertically into hyb tube using small hole, replace cap
- Incubate at 65°C overnight. Put wash solutions at 65°C overnight. Sign out of workstation and document isotope usage

# DAY 2



- ▶ **Washes:** Sign into workstation. Pour out hyb solution down sink using small hole in hyb tube (with lots of water). Wash blot sequentially, pouring old wash solutions down sink as above:
  - ▶ 50ml [Wash 1](#) for 10 mins
  - ▶ 50ml Wash 1 for 30 mins
  - ▶ 50ml [Wash 2](#) for 30 mins
    - ▶ For mismatched probes, adjust [Wash 2](#) for reduced [probe annealing temperature](#)
  - ▶ 50ml Wash 2 for 30 mins
    - ▶ While last wash is incubating, blank screen(s) with white light
- ▶ Blot off excess liquid using Whatman paper – do not allow membranes to dry completely  
Wrap in cling film and insert into cassette with screen, exposing to white side
- ▶ Leave for >24hours (depending on age of radiation) to expose
- ▶ Wash up all hyb tubes and leave to drain and sign out of workstation
- ▶ **Tidy up:** Dispose of tips, wash thoroughly with water/Decon and put in tub in bin
- ▶ Record radiation use/disposal on Isostock



# DAY 3

### ▸ Scan image:

- [Typhoon phosphorimager](#) in D69. Door code C148YZ
- Stick screen to solid metal tray white side up. Slide tray into phosphorimager (screen on underside).
- Open Typhoon software. Drag area to scan to correct size (usually J9).
- Pixel size 100  $\mu$ M. Sensitivity 1000 for colony hyb and 4000 for Southern.
- Choose name and location to save and start scan.
- Fill in Typhoon user book and radiation monitoring book (monitor bench with Geiger counter).
- Blank screen on light box
- Open image in Affinity photo to adjust levels

## PRE-HYB SOLUTION

Note: make solutions fresh in 50ml falcon

- ▶ 25ml dH<sub>2</sub>O
- ▶ microwave on high 10sec/tube
- ▶ 12ml 20X SSPE
- ▶ 2ml 20% SDS  
- ▶ 2ml 100X Denhardt's Solution

Total = 40ml

Pre-warm at 65°C in hyb oven

## 20X SSPE

- ▶ 175.3g NaCl
- ▶ 27.6g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$
- ▶ 7.4g EDTA (solid) or 40ml 0.5M EDTA (liquid)
- ▶ 800ml  $\text{dH}_2\text{O}$
- ▶ Adjust pH to 7.4 with NaOH
- ▶ make up to 1L with  $\text{dH}_2\text{O}$
- ▶ Autoclave





## 100X DENHARDT'S SOLUTION

- ▶ 1g Ficoll 400
- ▶ 1g PVP 360
- ▶ 1g BSA (fraction V)
- ▶ Make up to 50ml with dH<sub>2</sub>O
- ▶ Leave shaking in water bath until all solids are in solution
- ▶ Syringe filter 0.8µm
- ▶ Aliquot and store at -20°C

## HYB SOLUTION

Note: make solutions fresh in 50ml falcon

- ▶ 1.5g Dextran Sulphate (omit for colony hybs)
- ▶ ~18ml dH<sub>2</sub>O
- ▶ microwave on high 10sec/tube
- ▶ 9ml 20X SSPE
- ▶ 1.5ml 20% SDS



Total = 30ml

Pre-warm at 65°C in hyb oven

### WASH 1

- ▶ 50ml 20X SSPE
- ▶ 5ml 20% SDS
- ▶ dH<sub>2</sub>O up to 500ml



### WASH 2

- ▶ 5ml 20X SSPE
- ▶ 5ml 20% SDS
- ▶ dH<sub>2</sub>O up to 500ml



For mismatched probes,  
adjust Wash 2 for reduced  
[annealing temperature](#)

# PROBE ANNEALING TEMPERATURE VS SSPE/SDS CONCENTRATION

Length	% GC	%	[Na+]	%	Tm (°C)
100	65	100	50	0	80.3
25	65	100	50	0	62.3
50	65	100	50	0	74.3
75	65	100	50	0	78.3
100	65	100	50	0	80.3
125	65	100	50	0	81.5
150	65	100	50	0	82.3

X SSPE	% SDS	[Na+]	% GC	Length Tm (100 bp)	100 Tm (°C)
6	0.5	1009	65	98.4	98.4
2	0.5	348	65	93.0	93.0
0.5	0.5	100	65	85.1	85.1
0.5	0.1	86	65	84.0	84.0
0.2	0.5	50	65	80.3	80.3
0.2	0.1	37	65	78.2	78.2
0.2	0	33	65	77.4	77.4
0.1	0.5	34	65	77.6	77.6
0.1	0.1	20	65	73.8	73.8
0.1	0	17	65	72.7	72.7