

# RECOMBINATION ASSAY

# RECOMBINATION ASSAY – PRINCIPLE

- The assay is based on scoring recombination between a mutant chromosomal *leuB-Ag1* allele in a Haloferax strain (derivatives of H195 or H164 strains) and a mutant *leuB-Aa2* allele present on non-replicative pTA163 plasmid. Recombination between these alleles results in a wild-type *leuB+*, which is scored by the ability to grow on Hv-Min media lacking leucine.

## HALOFERAX TRANSFORMATION


### Day 1

- Use the general protocol for Haloferax transformation.

For each strain being assayed for recombination (including a WT control), two transformations are needed:

- Transform with 1 µg of **non-replicative** plasmid pTA168 (*pyrE2+*, *leuB-Aa2*) (dam- version of pTA163).
  - Plate transformants on **Hv-Min+Ura** (+Trp, +Thy, if relevant) at  **$10^{-1}$  -  $10^{-2}$**  dilutions to select for cells that have undergone recombination between the plasmid *leuB-Aa2* allele and the chromosomal *leuB-Ag1* allele, generating a wild-type *leuB+* allele.
  - To determine the total viable count, plate transformants on non-selective media (e.g., **Hv-YPC** or Hv-Ca-Ura (+Trp, +thy, if relevant) at  **$10^{-4}$  -  $10^{-6}$**  dilutions.
- Transform with 1 µg of **replicative** plasmid pTA357 (*pyrE2+*) (dam- version of p354).
  - Plate on Hv-Ca plates (+Trp, +Thy, if relevant) at  **$10^{-1}$  -  $10^{-2}$**  dilutions to select for cells that had taken up the replicative plasmid and thereby determine the transformation efficiency of the strain.
  - Plate on **Hv-YPC** at  **$10^{-4}$  -  $10^{-6}$**  dilutions to determine the total viable count.

The recombination frequency will be normalised to the transformation efficiency for each strain.



# RECOMBINATION ASSAY CONTINUED

## Day 5.

- Count colonies on plates and calculate recombination frequency (*RF*). This should be normalised to transformation efficiency (*TE*):
  - $RF = (Leu^{+}_{pTA168} / viable_{pTA168}) \times TE$ 
    - Where  $TE = (Ura^{+}_{pTA357} / viable_{pTA357})^{WT} / (Ura^{+}_{pTA357} / viable_{pTA357})^{Test}$
- Normalise recombination frequencies to WT, i.e.  $RF^{Test}$  relative to  $RF^{WT}$  (e.g. above, values in **bold**).
- Determine fraction of crossover (CO) and non-crossover (NCO) recombination events by patching  $Leu^{+}_{pTA168}$  colonies (from Hv-Min+Ura plates obtained in 1A) on:
  - Hv-Min** (+Trp, +Thy if relevant) to select for CO recombination events in cells that integrate pTA168, becoming *pyrE2+ leu+*
  - and then on **Hv-Min+Ura** (+Trp, +Thy if relevant), to ensure that all colonies patched are *leu+*
- CO events are *leuB+ pyrE2+*, the remaining recombination events (*leuB+, pyrE2-*) are NCO

Strain	WT	$\Delta hel308$	<i>hel308-F316A</i>	<i>hel308-D145N</i>	<i>hel308-D145N-F316A</i>
Recombination Frequency (RF)	1.92x10 <sup>-5</sup> (+/- 6.71x10 <sup>-6</sup> )	4.13x10 <sup>-4</sup> (+/- 2.73x10 <sup>-4</sup> )	1.81x10 <sup>1</sup> (+/- 4.59x10 <sup>0</sup> )	3.07x10 <sup>-5</sup> (+/- 2.05x10 <sup>-5</sup> )	9.35x10 <sup>-4</sup> (+/- 3.36x10 <sup>-4</sup> )
Transformation Efficiency (TE)	1.34x10 <sup>-3</sup>	4.85x10 <sup>-3</sup>	7.63x10 <sup>-3</sup>	4.23x10 <sup>-3</sup>	1.28x10 <sup>-1</sup>
Relative RF normalised by TE	1.43x10 <sup>-2</sup>	7.85x10 <sup>-2</sup>	2.37x10 <sup>3</sup>	7.26x10 <sup>-3</sup>	7.32x10 <sup>-3</sup>
	<b>1</b>	<b>5.5x</b>	<b>166000x</b>	<b>0.51x</b>	<b>0.51x</b>
CO fraction	16.5%	21.9%	<0.5%	35.7%	36.8%
NCO fraction	83.5%	78.1%	>99.5%	64.3%	63.2%

# PRACTICAL TIPS

- Over dry plates to avoid lawning. Dry plates for ~40 minutes as opposed to the standard 20 minutes.
- Always use freshly streaked cells from the -80°C and use within 1 week.